

Experiment #: 01

Experiment Title: 4 Point Method / Measurement of low resistances / Ohm's Law

Objectives:

1. To understand measure low resistances of materials
2. To correlate the method of **Four probe method** with low resistances of materials.

Four probe method is used to measure low resistances of materials.

Theory:

It is a **technique** to measure electrical impedance. It uses separate pairs for current carrying and voltage sensing electrodes, as shown in Fig.1.




Fig. 1 Four probe method

At a constant temperature, the resistance, R of a conductor is proportional to its length L and inversely proportional to its area of cross section A .

$$R = \rho \frac{L}{A} \quad (1)$$

Where ρ is the resistivity of the conductor and its unit is ohmmeter.

A semiconductor has electrical conductivity intermediate in magnitude between that of a conductor and insulator. Semiconductor differs from metals in their characteristic property of decreasing electrical resistivity with increasing temperature.

According to band theory, the energy levels of semiconductors can be grouped into two bands, valence band and the conduction band. In the presence of an external electric field it is electrons in the valence band that can move freely, thereby responsible for the electrical conductivity of semiconductors. In case of intrinsic semiconductors, the Fermi level lies in between the conduction band minimum and valence band maximum. Since conduction band lies above the Fermi level at $0K$, when no thermal excitations are available, the conduction band remains unoccupied. So conduction is not possible at $0K$, and resistance is infinite. As

temperature increases, the occupancy of conduction band goes up, thereby resulting in decrease of electrical resistivity of semiconductor.

Resistivity of semiconductor by four probe method

1. The resistivity of material is uniform in the area of measurement.
2. If there is a minority carrier injection into the semiconductor by the current- carrying electrodes most of the carriers recombine near electrodes so that their effect on conductivity is negligible.
3. The surface on which the probes rest is flat with no surface leakage.
4. The four probes used for resistivity measurement contact surface at points that lie in a straight line.
5. The diameter of the contact between metallic probes and the semiconductor should be small compared to the distance between the probes.
6. The boundary between the current carrying electrodes and the bulk material is hemispherical and small in diameter.
7. The surface of semiconductor material may be either conducting and non-conducting. A conducting boundary is one on which material of much lower resistivity than semiconductor has been plated. A non-conducting boundary is produced when the surface of the semiconductor is in contact with insulator.

Temperature dependence of resistivity of semiconductor

Total electrical conductivity of a semiconductor is the sum of the conductivities of the valence band and conduction band carriers. Resistivity is the reciprocal of conductivity and its temperature dependence is given by

$$\rho = A \exp\left(\frac{E_g}{2kT}\right)$$

Where E_g – band gap of the material

T – Temperature in kelvin

K – Boltzmann constant, $K = 8.6 \times 10^{-5} \text{ eV/K}$

The resistivity of a semiconductor rises exponentially on decreasing the temperature.

Applications


1. Remote sensing areas
2. Resistance thermometers
3. Induction hardening process
4. Accurate geometry factor estimation
5. Characterization of fuel cells bipolar plates

List of Equipment:

Apparatus


1. Measure voltage (V) **6**
2. supply current (A) to the surface of the crystal.
- 3.

The figure shows the arrangements of four probes that measure voltage (V) and supply current (A) to the surface of the crystal.



FYTONINX FOUR PROBE 2030

Circuit Diagram:



Procedure:

- Make the connections as per the circuit diagram.
- Connect the resistance wire in the left gap (between c & d) and resistance box in the right gap.
- Introduce some resistance in the circuit by taking out some resistance from the resistance box.
- Plug the key. Bring the jockey in contact with the end A first, and then with C. Note the deflection on the galvanometer.

Data Collection:

Table1:

W / S	f(W / S)
0.100	13.863
0.141	9.704
0.200	6.931
0.333	4.159
0.500	2.780
1.000	1.504
1.414	1.223
2.000	1.094
3.333	1.0228
5.000	1.0070
10.000	1.00045

Calculation:

Fig: 2 show the resistivity probes on a die of material. If the side boundaries are adequately far from the probes, the die may be considered to be identical to a slice. For this case of a slice of thickness w and the resistivity is computed as

$$\rho = \frac{\rho_0}{f\left(\frac{w}{s}\right)} \quad (2)$$

The function, $f(w/s)$ is a divisor for computing resistivity which depends on the value of w and s

We assume that the size of the metal tip is infinitesimal and sample thickness is greater than the distance between the probes,

$$\rho = \frac{V}{I} 2\pi s \quad (3)$$

Where V – the potential difference between inner probes in volts.

I – Current through the outer pair of probes in ampere.

S – Spacing between the probes in meter.

Result: