

## Experiment #: 08

### Experiment Title: Fresnel's biprism

#### Objectives:

1. To determine of the wavelength of light by interference with Fresnel's biprism.
2. Fresnel biprism is used to divide the wavefront of a monochromatic, coherent beam of light producing an interference pattern. The wavelength of the light is determined.

#### Theory:




Fig 1. Schematic representation of the beam path at Fresnel's biprism:

$A$ : Light source (He-Ne laser).

$A_1$  and  $A_2$ : Virtual light sources.

$S$ : Screen/wall.

$a$ : Distance between the two virtual light sources.

$d$ : Distance between two neighbouring intensity maxima or minima.  $P_1$  and  $P_2$ : Prism halves.

$L$ : Separation between laser and screen/wall.




Fig 2. Schematic representation of the image of the two virtual light sources:

H: Imaging lens H.

$A_1$  and  $A_2$ : Virtual light sources.

S: Screen/wall.

a: Distance between the two virtual light sources.

B: Distance between the image of two virtual light sources on the screen. g: Distance between the virtual sources and the lens H (object distance). b: Distance between the lens H and the screen S.

#### Precautions

1. Laser light is dangerous and can potentially cause visual impairment. Never look directly into any laser beam. Prolonged exposure will cause flash-blindness, after images and glare which will reduce or cause complete loss of visibility in the central field of vision.
2. Optical systems are sensitive and are often fine-tuned. Be very careful with the equipment, as a slight nudge might damage the equipment.
3. Stray light can obscure the images seen on the screen. Perform the experiment in pitch-black darkness.

#### List of Equipment:

#### Apparatus



- Fresnel biprism.
- Lenses ( $f = 20 \text{ mm}$  and  $300 \text{ mm}$ ). • Lens mounts.
- Swinging arm and slide mounts.
- Optical bench.
- Laser, He-Ne 1.0 mW, 220 V AC. • Measuring tape, 200 cm.

### Circuit Diagram:

### Procedure:

Part 1: Interference Pattern (Fig. 1):

In front of you, the He-Ne laser is mounted at the 2 cm mark.

Mount the lens ( $f = 20 \text{ mm}$ ) at the 23.3 cm mark. [This lens spreads the laser beam slightly \(widens it\).](#)

Mount the biprism at the 45 cm mark, with its tip facing the laser. [Use your finger to determine where the biprism's tip is.](#)

You should be able to see an interference pattern on the wall.

Using a vernier scale, measure the separation  $D$  between five maxima/minima. Repeat step 5 three times, measuring  $D_1$ ,  $D_2$  and  $D_3$ .

Measure the distance  $L$  between the laser and the screen/wall.

Part 2: Virtual Source Separation (Fig. 2):

Mount the lens ( $f = 300$  mm) at  $\approx 60$  cm mark. You should be able to see two separate light points.

Using a vernier scale, find the distance  $B$  between the two light points. Measure the distance  $b$  between the lens H and the screen S.

7 Evaluation

1. Calculate the separation  $d$  between two successive maxima/minima:

$$d_n = \frac{D_n}{5}$$

Where  $n = 1, 2$  and  $3$ .

2. Find the average  $d_{avg}$  of  $d_n$ .

3. Calculate the distance  $g$  between the virtual sources and the lens H (object distance) using the imaging equation:

$$g = \frac{fb}{b-f}$$

Where  $f = 300$  mm.

Calculate the distance  $a$  between the virtual light sources:

$$a = \frac{Bg}{b}$$

Find the wavelength  $\lambda$  using:

$$\lambda = \frac{d_{avg}a}{L}$$

Calculate the error percentage.

Postlab Questions

What is meant by the virtual source in Fresnel's Biprism experiment? Using the experimental sketch above, explain how such a source arises.

You can find the value of  $g$  using an equation other than the imaging equation. Write the expression of that equation, explaining how you obtained it.

**Data Collection:**

**Calculation:**

**Result:**