FYTRONIX
R&D GROUP TECHNOLOGIES

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


Cyclic voltammogram of a 40 nm TiO_2 electrode acquired between 0.8 V and 3.2 V at scan rates between 0.5 mV/s and 5 mV/s.

Potential evolution during constant-current charging and discharging of a 40 nm TiO_2 electrode at a current density of ± 0.18 µA/cm^2.


EIS for Supercapacitor Electrodes

Electrochemical Impedance Spectroscopy (EIS) is a powerful electrochemical technique used to analyze resistance, charge transfer processes, and ion diffusion in supercapacitor electrodes.


From Nyquist:

- X-axis: Real impedance (Z')
- Y-axis: Imaginary impedance (-Z")

This plot provides information about the electrochemical processes occurring at different frequencies.

Main Regions in the Nyquist Plot


1. High-Frequency Region


The intercept of the curve with the real axis represents the solution resistance (Rs). Rs includes contributions from:

- Electrolyte resistance
- Electrode intrinsic resistance
- Contact resistance between electrode and current collector

Key meaning:

- Lower Rs → better electrical conductivity
- Higher Rs → larger internal resistance in the system

Imaginary Impedance (-Z'')
Rs
(Solution
Resistance)
Semicircle
Arc (Rct)
Rct
WDR
Warburg
Impedance
(WDR)
High Frequency
Region
Low Frequency
Region

2. Semicircle Region (Charge Transfer Region)


A semicircle usually appears at high to mid frequencies.

- The diameter of the semicircle represents the charge transfer resistance (Rct).
 What the semicircle indicates:

- Large semicircle → high Rct → slow electron transfer
- Small semicircle → low Rct → fast electrochemical reactions
- Smaller semicircle → better electrode kinetics


Interpretation:

- Small semicircle = efficient charge transfer
- Large semicircle = poor electrochemical activity

3. Straight Line Region (Low Frequency Region)

At low frequencies, the Nyquist plot becomes a straight line. This line is related to Warburg diffusion resistance (WDR), which reflects ion diffusion within the porous electrode material.

Compensation
Start:
10 Hz
Steps:
201
Mode:
Amplitude
Stepping:
Logarithmic
Enable Open
RUN
Stop:
10 MHz
Offset:
0 V
Amplitude:
1 V
Resistor
10 kΩ
Enable Short
ABORT
|Z|
|Y| (S)
PHASE
CS (F)
CP (F)
D
LS (H)
LP (H)
Q
RS
G (S)
RP
X (OHMS)
B (S)
Vrms
Irms
RS - X
Compensation
Reset
Open
Save
Linear
1E+1
8E+0
6E+0
4E+0
2E+0
0E+0
| Z | (OHMS)
1,00
11,00
21,00
31,00
41,00
51,00
61,00
71,00
81,00
91,00
101,00
Freq (Hz)
Series1
Check Resistor Value!